


8D

3 IN 1 SYSTEM
H I F U


3 HANDLES AVAILABLE



Work Theory



8D HIFU is Non Surgical Face Lift Treatment- 8D HIFU (High Intensity Focused Ultrasound) Technology , the aging effects will show on the skin like the wrinkles and the eye-bags. It is the signals to alert you your collagen and elastin are losing from your body. The collagen is the key structural component to support your skin.



8D HIFU (High-Intensity Focused Ultrasound) which is the new technical for face lifting and tightening. 8D HIFU is delivered into the deep skin and cause thermal coagulation to SMAS (Superficial Muscular Aponeurotic System). It can provide the skin tightening and rejuvenation happens gradually with the time. 8D HIFU technology device gives effective, safe and painless skin tightening.

RF TIGHTENING



Strengthen and tighten relaxed skin with remarkable tightening effect

Only one treatment is needed to smooth wrinkles and improve them comprehensively.


Make skin texture more delicate and glossy

Filling in shrinkage and loss of collagen vacancies to enhance the relaxation profile

Remove the traces of age such as roughness, and Skin rejuvenation

Make skin tender, white and youthful

Specifications



Frequency

4MHz

Power

350W

Treatment Depth

DS1.5mm DS3.0mm DS4.5mm (Standard)
DS6.0mm DS8.0mm DS10.0mm
DS13.0mm DS16.0mm (Optional)

Treatment Length

8-25mm

Width

0--10mm

Energy

0.1-- 2.0J

Pitch

1.0mm-2.5mm

Clearanceh

1.0mm-2.0mm

Power Source

AC 110V~220V--50Hz/60Hz

System introduction

①

Single Line

②

Multi
Line

③

Single Dot

④

12.1 inch screen

⑤

Switch

⑥

Trolley



Advantages


- ▶ Completely non-invasive and non-ablative to the epidermis layer
- ▶ No down time, and make-up can be wearied only in 15 min after the treatment
- ▶ Fast operation time , full face and neck only takes 40-60 minutes
- ▶ Long lasting result, 2-3 year can last only with one treatment
- ▶ Flexible treatment, one single treatment achieves instant result. High quality, low consumables





Treatment: Body

Multi Line points: 10000



Treatment: Face, Neck, Jaw

Single Line points: 10000

2019



Treatment: Around eyes and mouth

Single Dot Points: 60000

Multi Line Options



1.5 mm



3.0 mm



4.5 mm



6.0 mm



8.0 mm



10.0 mm



13.0 mm



16.0 mm

Single Line Options



1.5 mm



2.0 mm



3.0 mm



4.5 mm



6.0 mm



9.0 mm



13.0 mm

Single Dot Options



1.5 mm



3.0 mm



4.5 mm



8.0 mm



13.0 mm



8.0 mm



13.0 mm



1.5 mm



3.0 mm



4.5 mm



1.5 mm



3.0 mm



4.5 mm



Q&A

Who is suitable for 8D HIFU treatment ?

It suit for 30 to 60 years old customer who want to reduce the pore size and looks younger.

when you will see the treatment result?

It suit for 30 to 60 years old customer who want to reduce the pore size and looks younger.

Q&A

why some clients have good effect of face lifting, but some clients have hollow cheek after treatment?

80% effect comes from the after treatment recovery, it's depend on patients skin property and recruitment of cell-repairing factor (such as DPHP)

(If patients age more than 30 years old, the collagen will loss more and more seriously, so before treatment, they need drink some Special repair of collagen ACMETEA so that can help regeneration of collagen fibers)

During these time, If Patients didn't support enough collagen timely, due to the deficiency of regeneration collagen will result in hollow cheek.

Q&A


If a patient undergo 8D HIFU treatment too frequently, what risk he will take ? what is the suitable interval time for 8D HIFU treatment?

Usually the treatment results of 8D HIFU treatment can last 2-3 years. 8D HIFU treatment has a little damage on Subcutaneous, so we don't suggest a patient accept 8D HIFU treatment too frequently, we suggest the interval time is above 2 years.


How to do the postoperative care after 8D HIFU treatment?

After the treatment, the patient can do the daily life. There is no special considerations. The treatment area maybe swelling, but don't worry, it will disappear after a few hours.


UI Interface




Single Line



Single Dot



Multi Line



Multi Line



BEFORE & AFTER




Before

After



Before

After



Before

After



Before

After



Before

After



Before

After